

Handy Guide

Basic Clinical Photography

An excellent clinical image can make all the difference in a telemedicine consultation. Here are some tips to help make your images the best that they can be.

SUBJECT: Should fill the majority of frame. Also keep area free of clutter that may distract from the image.





SCALE: Use a ruler or other scale to provide information about the size of the subject.



NOTE: Use solid colored background, but avoid white or black. Also avoid textures and reflective surfaces.











PERSPECTIVE: Keep camera parallel to the subject surface.





MACRO: Use for all close-ups (2"-28"). Avoid zoom when possible as it makes it harder to hold the camera steady.







AFHCAN AFHCAN www.afhcan.org

Handy Guide

Basic Clinical Photography

FOCUS: On most cameras, gently squeeze shutter to auto-focus, then depress shutter completely to take the picture.





STEADY THE CAMERA:

Use your arms as a tripod on a flat surface or against your sides.

Gently squeeze shutter.





FLASH: Try photos with and without flash to see which produces the best image.





GOAL

Subject clear



Focus sharp



Scale shown



Well lit



True color



Great Image!





The Core Standards for Telemedicine Operations as well as the practice guidelines for Teledermatology can be found on ATA's website: http://www.americantelemed.org